Nord FLOW

EN	Light curing nano flowable composite
DE	Lichthärtendes nano-fließfähiges Komposit
BG	Лекo втвърдяващ се нанотечен композит
FR	Composite nano-fluide photopolymérisable
EL	Ελαφρύ σκληρυνση nano -ρευστό σύνθετο
IT	Composito fluido fotopolimerizzabile nano
PL	Światłoutwardzalny nano płynny kompozyt
ES	Composite nano fluido fotopolimerizable
TR	İşıkla sertleşen nano akışkan kompozit

INSTRUCTION FOR USE EN

DESCRIPTION
Nord FLOW is a light curing nano flowable, bioinert, radiopaque composite under the Vita® shades. Nord FLOW is a high aesthetic, highly resistant, superior polishability product with optimal flow characteristics.

COMPOSITION
Dental glass grinded 50-70%, methacrylate mixture 30-40%, silicon dioxide 1-5%, coinitiator <1%, photoinitiator <1%, stabilizer <1%, inhibitor <1%, opacifier <1%, pigment 1%.
Nord FLOW does not contain medicinal substance, including human blood or plasma derivative; tissues or cells, or their derivatives, of human origin; or cells of animal origin, or their derivatives, as referred to in Regulation (EU) No. 722/2012; substances which are carcinogenic, mutagenic, toxic to reproduction or having endocrine-disrupting properties.

PERFORMANCE CHARACTERISTICS

light curing	20-30s
depth of cure	3.24±0.03 mm
flexural strength	110.2±4.1 MPa
water sorption	17.10±0.20 µg/mm ³
water solubility	0.00±0.00 µg/mm ³

INTENDED PURPOSE AND CLINICAL BENEFITS
Nord FLOW restores/improves aesthetic appearance of restorable tooth; restores/maintains dental function of restorable tooth; protects biological structures of restorable tooth and adjacent tissues.

CLINICAL INDICATIONS

- For restorations of class III, IV and V cavities; root surface caries restorations;
- For sealing pits and fissures;
- For initial placement in class I and II cavities.

CONTRA-INDICATIONS
Patients who have a history of severe allergic or irritation reactions to product or to any of the ingredients.

RESTRICTIONS TO COMBINATIONS
Nord FLOW should not be used with products containing eugenol as eugenol may disturb polymerization process.

UNDESIRABLE SIDE EFFECTS
In susceptible individuals, Nord FLOW may cause allergic or irritation reactions (skin, eye, mucosa, respiratory tract irritation).

RESIDUAL RISKS
Risk control measures have been implemented and verified, risk is reduced as far as possible, the overall residual risk is judged to be acceptable.

PATIENT TARGET GROUP
No restrictions known regarding patient population, their age and general health conditions. There may be children, middle aged or elderly patient.

INTENDED PART OF THE BODY OR TYPES OF TISSUES OF BODY FLUIDS
Part of the body – mouth. Tissues or body fluids contacted by the device – tooth, saliva.

INTENDED USER
Nord FLOW is developed for professional use in dentistry only. Its user only licensed doctor who has knowledge how to use common dental composites. There is no need for specific training.

STERILITY
Nord FLOW is delivered non-sterile. There is no need of any preparatory sterilization, cleaning or disinfection, preventive, regular maintenance or calibration to ensure that the device operates properly and safely during its intended lifetime. However, do not use if primary package is damaged.

USE ENVIRONMENT
Nord FLOW is designed to be used in dental office where ambient temperature is 18-25°C. Dispensed amount of composite is suitable for single use (only for one patient). Do not re-use. Dispensed amount kept not in original package may lead to loss of function.

CONSUMABLE COMPONENTS AND ACCESSORIES
No accessories are supplied with the device. Consumables, such as application tips, are supplied with the device.

INSTRUCTION FOR USE

CAVITY PREPARATION:

1. Prepare cavity as always.
2. Clean the surface with oil-free prophylaxis paste, such as Nord FASTE.
3. For deep cavities use calcium hydroxide liner or glass ionomer base lining cement.

ETCHING, BONDING:

1. Apply layer of etch, such as Nord GEL to surface to be etched. Leave etch in place for 15 seconds (dentine), 30 seconds (enamel). Rinse with water and dry with air. Avoid over drying dentin.
2. Apply a layer of adhesive, such as Nord BONDING LC immediately onto etched surface, follow manufacturer's instruction for use.

SYRINGE PREPARATION:

1. Remove syringe cap.
2. Promptly and carefully attach the dispensing tip to the syringe.
3. Test flow of materials from tip before using intraorally.

PLACEMENT OF Nord FLOW:

1. Before bringing the syringe to the mouth, remove the air from the dispensing tip. To remove air from the tip, with the tip pointing upwards, gently push forward the syringe plunger. If the air is still inside the dispensing tip, air bubbles may be removed at the time of injection.
2. Delicate push on plunger and apply layer of material into the cavity. Do not force plunger.
3. Do not apply layers more than 2 mm deep.
4. Light cure for 20-30 seconds (depends on layer deep). Use LED polymerization lamp with light intensity 1200mW/cm² in full mode, not ramp or pulse mode. Some lamps with higher intensity could require less time of polymerization, follow manufacturer's instruction for use. Finish restoration.

WARNINGS
After the desired amount of material extruded, immediately remove application tip and close the syringe cap, so that the material is not unlighted. The material is sensitive to light. Avoid too long manipulation time under intensive lighting. Do not use Nord FLOW for patients who have a history of severe allergic or irritation reactions to product or to any of the ingredients. Nord FLOW does not emit radiation and does not cause any electromagnetic interferences.

PRECAUTIONS
It is recommended to use cofferdam during application of the product. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists: Get medical advice/attention. IF ON SKIN OR MUCOSA: Wash with plenty of water. If skin/mucosa irritation or rash occurs: Get medical advice/attention. Take off contaminated clothing and wash before reuse. IF SWALLOWED: Rinse mouth. Call a Poison Center or doctor/physician if you feel unwell. IF INHALED: Remove person to fresh air and keep comfortable for breathing. Wash hands thoroughly after handling. Use only in a well-ventilated area. It is recommended to wear protective gloves/protective clothing/eye protection/face protection for doctor and patient.

PRECAUTIONS to be taken in the event of changes in the performance of the device:
If noticed any changes in the performance of the device, immediately remove applied material with suitable dental instrument. Ask patient how she/he is feeling. If patient noticed any undesirable side-effects, immediately call to a local poison center.

SHELF-LIFE
Shelf-life of Nord FLOW is 4 years from the date of manufacture. Do not use after the expiry date. The batch number should be quoted in all correspondence. See packaging for batch and expiry date.

STORAGE
Keep product tightly closed in dry well-ventilated place at 4-28°C. Protect from direct sunlight and heat sources. Do not freeze. Keep out of the reach of children!

DISPOSAL
Dispose of contents/container to as required by national regulatory requirements.

VIGILANCE
If any serious incident that has occurred in relation to the device report to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

SUMMARY OF SAFETY AND CLINICAL PERFORMANCE
Nord FLOW is safe and performs as intended if it is used in accordance to manufacturer's instruction for use. Summary of safety and clinical performance will be introduced in EUDAMED as soon as it will start work.

MANUFACTURERS RESPONSIBILITY
Our products have been developed for professional use in dentistry. As the application of our products is beyond our control, the user is fully responsible for the application. Of course, we guarantee the quality of our products in accordance with the applied standards.

VALIDITY
Upon publication of this instruction for use, all previous versions are superseded.

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* Registered trademark of the Vita Zahnfabrik H.Rauter GmbH & Co. KG, Bad Sackingen, Germany.
* Vita Zahnfabrik H.Rauter GmbH & Co.nun (KG, Bad Sackingen, Almanya) tescilli ticari markasıdır.