Nord PAK

EN	Knitted retraction cord
DE	Gestrickte Retraktionsschnur
BG	Зъболекарски конец
FR	Cordon de rétraction tricoté
EL	Πλεγμένο νήμα ανώθησης
IT	Cordone di retrazione a maglia
PL	Dzianinowa nić retrakcyjna
ES	Cordón de retracción de punto
TR	Dokuma retraksiyon kordonu

INSTRUCTION FOR USE EN

DESCRIPTION

Nord PAK is non-sterile cord made of 100% cotton, knitted into thousands of tiny loops forming long interlocking chains. Nord PAK is available in 4 sizes: #1, #0, #00, #000. Marking depends on the thickness of cord. A unique attribute of this cord is that it doesn't entangle in the diamond bur.

COMPOSITION

Cotton 100%.
 Product does not contain medicinal substance, including human blood or plasma derivative; tissues or cells, or their derivatives, of human origin; tissues or cells of animal origin, or their derivatives, as referred to in Regulation (EU) No. 722/2012; substances which are carcinogenic, mutagenic, toxic to reproduction or having endocrine-disrupting properties.

PERFORMANCE CHARACTERISTICS

Colour	knitted retraction cord 000 - black knitted retraction cord 00 - yellow knitted retraction cord 0 - lilac knitted retraction cord 1 - blue
Diameter	knitted retraction cord 000 - 1.0mm (ultra-thin) knitted retraction cord 00 - 1.2mm (thin) knitted retraction cord 0 - 1.35mm (medium) knitted retraction cord 1 - 1.5mm (medium thick)

INTENDED PURPOSE AND CLINICAL BENEFITS

Nord PAK is intended for use during dental procedures for gingival retraction.

CLINICAL INDICATIONS

- **#1 for the anterior teeth preparation;**
- **#0 for anterior teeth; as upper cord for "two cord" technique;**
- **#00 during preparation and cementation of veneers; restorative procedures dealing with thin, friable tissues;**
- **#000 for lower anteriors; when luting near gingival and subgingival veneers; class III, IV V restorations; first cord for "two cord" technique.**

CONTRA-INDICATIONS

Patients who have a history of severe allergic or irritation reactions to product or to any of the ingredients.

RESTRICTIONS TO COMBINATIONS

None known.

UNDESIRABLE SIDE EFFECTS

In susceptible individuals, product may cause allergic or irritation reactions (skin, eye, mucosa, respiratory tract irritation). Do not use product if patient is allergic to product or to any of the ingredients.

RESIDUAL RISKS

Risk control measures have been implemented and verified, risk is reduced as far as possible, the overall residual risk is judged to be acceptable.

PATIENT TARGET GROUP

No restrictions known regarding patient population, their age and general health conditions. There may be children, middle aged or elderly patient.

INTENDED PART OF THE BODY OR TYPES OF TISSUES OF BODY FLUIDS

Part of the body – mouth. Tissues or body fluids contacted by the device – tooth, saliva, mucosa.

INTENDED USER

Nord PAK is developed for professional use in dentistry only. Its user only licensed doctor who has knowledge how to use common dental materials. There is no need for specific training.

STERILITY

Nord PAK is delivered non-sterile. There is no need of any preparatory sterilization, cleaning or disinfection, preventive, regular maintenance or calibration to ensure that the device operates properly and safely during its intended lifetime. However, do not use if primary package is damaged.

USE ENVIRONMENT

Nord PAK IS designed to be used in dental office where ambient temperature is 18-25°C. Do not re-use.

CONSUMABLE COMPONENTS AND ACCESSORIES

No consumable components and accessories are supplied with the device.

INSTRUCTION FOR USE

1. Prepare and dry the tooth.

2. Choose the cord size. Because cord compresses upon packing, use a cord size that appears too large.
3. Cut the required length of retraction cord.
4. Soak the cord in hemostatic liquid/gel (not applicable for impregnated cord).
5. After preparing and drying the tooth, retraction cord is placed closely around the cervix and inserted carefully into the sulcus using a small retraction cord packer or spatula. The thin cord packer quickly slips it into position.
6. During the following preparation for the procedures, the mechanical effect of the cord and chemical effect of hemostatic material, results in a definite retraction of the gingiva.
7. Perform the necessary procedures.
8. If any bleeding is present, repeat steps until hemostasis is achieved.
9. When the procedures are completed, remove cord and rinse with water.

WARNINGS

Do not use product for patients who have a history of severe allergic or irritation reactions to product or to any of the ingredients. Product does not emit radiation and does not cause any electromagnetic interferences.

PRECAUTIONS

Wear protective gloves/protective clothing/eye protection/face protection for doctor and patient.

SHELF-LIFE

Shelf-life of Nord PAK is 6 years from the date of manufacture. Do not use after the expiry date. See packing for expiry date. Do not re-use.

STORAGE

Keep product tightly closed in dry well-ventilated place at 4-28°C. Protect from direct sunlight and heat sources. Do not freeze. Keep out of the reach of children!

DISPOSAL

Dispose of contents/container to as required by national regulatory requirements.

VIGILANCE

If any serious incident that has occurred in relation to the device report to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

SUMMARY OF SAFETY AND CLINICAL PERFORMANCE

Nord PAK is safe and performs as intended if it is used in accordance to manufacturer's instruction for use.

MANUFACTURERS RESPONSIBILITY

Our products have been developed for professional use in dentistry. As the application of our products is beyond our control, the user is fully responsible for the application. Of course, we guarantee the quality of our products in accordance with the applied standards.

VALIDITY

Upon publication of this instruction for use, all previous versions are superseded.

PACKAGING

- REF NRDPU11** Knitted retraction cord #1 (244cm), medium thick, blue
- REF NRDPU01** Knitted retraction cord #0 (244cm), medium, lilac
- REF NRDPU02** Knitted retraction cord #00 (244cm), thin, yellow
- REF NRDPU03** Knitted retraction cord #000 (244cm), ultra-thin, black

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